












Variolink® Esthetic Adhesive Cementation System

Shade reference table Variolink® Esthetic | Variolink® II

Variolink II						
Shade designation	White opaque	Bleach XL-010	Transparent	White 110/A1	Yellow 210/A3	Brown 340/A4
Product Code Refill Base	558954	558949	558950	558951	558952	558953
Variolink Esthetic						
Shade designation	Light+	Light	Neutral	Warm	Warm+	
Product Code Refill LC	666126	666127	666128	666129	666130	
Product Code Refill DC	666117	666118	666119	666120	666121	

The shades of Variolink Esthetic slightly differ from those of Variolink II.

In clinical conditions, the Variolink Esthetic shades achieve the same shade effect underneath restorations as the corresponding Variolink II shades.



FRIENDSHIP DENTAL LABORATORIES, INC.

8510 Philadelphia Road
Rosedale, MD 21237

Phone: 410-780-7700













Email: info@friendshipdentallab.com

WWW.FRIENDSHIPLABS.COM



Variolink® Esthetic Adhesive Cementation System

Shade reference table Variolink® Esthetic | Variolink® Veneer

Variolink Veneer							
Shade designation	High Value +3	High Value +2	High Value +1	Medium Value 0	Low Value -1	Low Value -2	Low Value -3
Product Code Refill Base	591423	591424	591425	591426	591427	591428	591429
Variolink Esthetic							
Shade designation	Light+	Light	Neutral	Warm	Warm+		
Product Code Refill LC	666126	666127	666128	666129	666130		

The shades of Variolink Esthetic slightly differ from those of Variolink Veneer. In clinical conditions, the Variolink Esthetic shades achieve the same shade effect underneath restorations as the corresponding Variolink Veneer shades.



FRIENDSHIP DENTAL LABORATORIES, INC.

8510 Philadelphia Road

Rosedale, MD 21237

Phone: 410-780-7700

Email: info@friendshipdentallab.com

WWW.FRIENDSHIPLABS.COM

